

FORCEPS DECON STATION

Another Trusted ViraTek Product




How do you decontaminate your forceps when changing cages?

- Do you place them in an open container of chemical?
- Did you know that chemicals such as chlorine dioxide “gas off” in less than 2 hours rendering them ineffective?

Do you give the forceps enough time to be decontaminated?

- Did you know that most decontaminating chemicals require a 10 minute contact time?
- Did you know that if you just dip forceps in a container of chemical, they will not receive the necessary contact time

The new FORCEPS DECON STATION by ViraTek, Inc. addresses all of these issues



- Easily fits under a hood
- Maintains the efficacy of the disinfectant for at least 8 hours
- Reduces the usage of disinfectant
- Increases disinfecting contact time
- Eliminates spills
- Ensures forceps rotation through the use of numbered holders
- Durable PVC construction stands up to chemicals

Contact ViraTek for a quote



ViraTek, Inc.

www.viratek.com • info@viratekinc.com

13617 Silver Hill Road
Sumerduck, VA 22742
Ph 540-905-5284
Fax 540-439-4048